

Benchtop Centrifuge

2-5

OPERATING MANUAL

Dear customer,

Congratulations for purchasing a SIGMA laboratory centrifuge. You have selected a device which combines many advantages.

The electronic operation control allows a trouble-free use of the centrifuge. With its 3-phase drive, maintenance-free quiet operation without any carbon dust pollution is guaranteed.

Your device is equipped with user-friendly options which make the operation and standard settings easier for you. Built-in error-detecting functions keep the user from entering incorrect values and check the complete operation.

All settings are executed via the control panel with a coated surface which protects the device against moisture and dust. In addition, the interior of the centrifuge is also easy to clean. We offer you a device that combines functional variety with practical applications.

We thank you for your confidence and wish you a successful application of the centrifuge.

SIGMA Laborzentrifugen GmbH
Postfach 1713 - D-37507 Osterode
Tel. 05522/5007-0 - Telefax 05522/500712
Internet: www.sigma-zentrifugen.de
eMail: info@sigma-zentrifugen.de
SIGMA Service-Tel. 05522/5007-25

Konformitätserklärung
(73/23/EWG; 89/336/EWG; 98/37/EWG)
Statement of Conformity
(73/23/CEE; 89/336/CEE; 98/37/CEE)
Déclaration de conformité
(73/23/CEE; 89/336/CEE; 98/37/CEE)

Die nachfolgend bezeichnete Maschine wurde in Übereinstimmung mit den Richtlinien 73/23/EWG; 89/336/EWG und 98/37/EWG hergestellt und geprüft.

The following machine is manufactured and tested in compliance with directions 73/23/CEE; 89/336/CEE and 98/37/CEE.

La machine désignée ci-dessous est produit et examiné conforme aux directives 73/23/CEE; 89/336/CEE et 98/37/CEE

Bezeichnung der Maschine: Laborzentrifuge
Machine: Laboratory Centrifuge
Désignation de la machine: Centrifugeuse de laboratoire

Maschinentyp : 2 - 5
Type:
Type de la machine:

Bestell Nr. : 10134, 10135, 10136
Part No.:
Réf. usine:

Normen: EN 61010-2-020
Standards: EN 61000-3-2 ; EN 61000-3-3
Normes : EN 61326

Sigma Laborzentrifugen
An der Unteren Söse 50
D-37520 Osterode



01.02.2002

Geschäftsführer
Managing Director
Directeur Gérant

.....
Fabr. Nr. Serial No. Numéro de fabrication

General Information

- 1.1 Technical data
- 1.2 Suitable accessories
- 1.3 Scope of supply
- 1.4 Standards and regulations
- 1.5 Safety instructions
- 1.6 Symbol table

Description of the Centrifuge

- 2.1 General outlay
- 2.2 Construction and constructive safety measures
- 2.3 Drive
- 2.4 Operation and display
- 2.5 Electronics
- 2.6 Safety devices
 - 2.6.1 Lid lock, cover closing device
 - 2.6.2 Standstill monitoring
 - 2.6.3 System check
 - 2.6.4 Ground wire check
 - 2.6.5 **Option:** Imbalance monitoring system

Installation and Start-up

- 3.1 Unpacking of the centrifuge
 - 3.1.1 Transport safety device
- 3.2 Installation
 - 3.2.1 Site
 - 3.2.2 Connection
 - 3.2.3 Fuses / emergency circuit breaker on site
- 3.3 Installation of rotors and accessories
- 3.4 Initial start-up
 - 3.4.1 Switching on of the centrifuge
 - 3.4.2 Opening lid
 - 3.4.3 Installation of an angle rotor and a swing-out rotor

Operating Elements

- 4.1 Operating panel
 - 4.1.1 Start-key
 - 4.1.2 Stop-key with soft and softstop function
 - 4.1.3 Lid-key
 - 4.1.4 Knobs

- 4.2 Displays
- 4.2.1 Speed/RCF (Relative centrifugal force)
- 4.2.1.1 Speed
- 4.2.1.2 Relative centrifugal force (RCF)
- 4.2.2 Time
- 4.2.3 Standard run mode
- 4.2.3.1 Continuous run mode
- 4.2.3.2 Short-time operation
- 4.2.3.2.1 Short-run 1
- 4.2.3.2.2 Short-run 2

Notes for Centrifugation

- 5.1 Practical notes for centrifugation
- 5.2 Forbidden centrifuging operations

Care and Maintenance

- 6.1 Care and cleaning of centrifuge
- 6.2 Care and cleaning of accessories
- 6.2.1 Service notes
- 6.3 Rotor trunnion pins
- 6.4 Glass breakage
- 6.5 Sterilization and disinfection of rotor chamber and accessories
- 6.5.1 Autoclaving
- 6.6 Checks by operator

Appendix

- 7.1 Mathematical relations
- 7.1.1 Relative Centrifugal Force (RCF)
- 7.1.2 Density
- 7.2 Error correction
- 7.2.1 Centrifuge cannot be started
- 7.2.2 Centrifuge decelerates during operation
- 7.2.3 Lid cannot be opened/closed
- 7.2.4 Emergency lid release
- 7.2.5 Problems with the centrifuge?
- 7.3 Error mode
- 7.3.1 Error codes
- 7.4 Diagram speed/gravitational field
- 7.5 Declaration of decontamination/Return declaration
- 7.6 Leaflet

1.1 Technical Data	
Manufacturer:	S I G M A Laborzentrifugen GmbH 37520 Osterode
Type:	2-5
Electr. connection: Protection class:	see nameplate I
Power consumption (kVA): Rated power (kW): Max. current (A):	0,17 0,13 0,8 (230 V/50 Hz) respectively 1,6 (120 V/60 Hz)
Power data:	
Max. speed (rpm): Max. capacity (ml): Max. gravitational field (x g): Max. kin. energy (Nm):	4 000 400 2 576 1 694
Further parameters	
Time range:	0 - 30 min/continuous run/ short-time operation
Dimensions:	
Depth (mm): Width (mm): Height (mm): Weight (kg): EMC (acc. to EN 55011): Noise level (dBA):	452 365 300 21 Class B < 58
Notes of user:	
Serial number:
Supply date:
Inventory number:
Location:
Responsibility:

The figures are valid for an ambient temperature of 23 °C +/- 2 °C and nominal voltage 220 – 240 volt +/- 5 %. (Allowable ambient temperature +4 °C - +40 °C; max. humidity 80 %.)
Subject to technical alterations.

1.2 Accessories Suitable for SIGMA 2-5

Part No.	Description	Max. speed (rpm)	Max. gravitational field (x g)
11031	Swing-out rotor 16 x 15 ml, complete, consisting of rotor 11030, 4 carriers 13020 and 16 polystyrene tubes 15020, max. radius 13.6 cm, min. radius 5 cm	4 000	2 433
11032	Swing-out rotor 4 x 100 ml, complete, consisting of rotor 11030, 4 buckets 13099 and 4 glass tubes 15100, max. radius 13.4 cm, min. radius 5 cm	4 000	2 397
11030	Swing-out rotor for 4 buckets or carriers 13020, 13031, 13033, 13034, 13035, 13038, 13099		
11035	Swing-out rotor 8 x 15 ml, complete, consisting of rotor 11036, 2 carriers 13020 and 8 polystyrene tubes 15020, for tubes max. Ø 17 x 100 mm, e.g. Monovettes, max. radius 13.4 cm, min. radius 3.7 cm	4 000	2 397
11036	Swing-out rotor for 2 buckets or carriers 13020, 13031, 13033, 13034, 13035, 13038, 13099		
13020	Multiple carrier for 4 round or conical tubes 15 ml, max. Ø 17 x 90 - 110 mm, e.g. 15015, 15020, 15022, 15023, 15024, Monovettes, suitable for 11030, 11032, 11036, polyallomer		
17907	Adapter for 1 glass tube 6 ml, e.g. 15007, 15027 or 1 Monovette 4.5 ml Ø 11 x 108 mm, suitable for 13020		
13035	Multiple carrier for 4 round or conical tubes 5-15 ml, max. Ø 17 x 80 - 110 mm, e.g. 15000, 15010, 15015, 15020, 15022, 15023, 15024, 15039, Monovettes, suitable for 11030, 11031, 11032, 11035, 11036, polyallomer		

Part No.	Description	Max. speed (rpm)	Max. gravitational field (x g)
13099	Bucket for 1 tube 100 ml, e.g. 15100, 15102, 15103, 15106, suitable for 11030, 11031, 11035, 11036, polyallomer		
17950	Adapter for 1 tube 50 ml, e.g. 15049, 15050, 15056, 15057, suitable for 13099		
17925	Adapter for 1 glass tube 25 ml, e.g. 15025, 15026, suitable for 13099		
13033	Multiple carrier, aluminium, for 5 round or flat bottom tubes approx. 5 ml, max. Ø 12.8/16.5 x 65 - 85 mm, e.g. 15060 (RIA), hemolyse or Vacutainer tubes, incl rubber cushion 16005, suitable for 11030, 11031, 11032, 11035, 11036		
13034	Bucket, aluminium, incl. polysulfone screw cap 17130, for round carriers 14029, 14030, 14031, 14032, 14033, 14034, 14035, max. tube length 110 mm, suitable for 11030 (from year of construction 2003), 11031, 11032, 11035, 11036, max. radius 14.1 cm, min. radius 7.1 cm	4 000	2 522
14029	Round carrier for 5 round bottom tubes 5-7 ml, max. Ø 12.5/15 x 65 - 105 mm, e.g. 15007, 15027, 15060, polypropylene, suitable for 13034		
14030	Round carrier for 4 round bottom tubes 10-12 ml, max. Ø 16.2/17.5 x 80 - 110 mm, e.g. 15000, 15010, 15015, 15024, 15039, polyethylene, suitable for 13034		
14031	Round carrier for 1 round bottom tube approx. 25 ml, max. Ø 25,3/30 x 70 - 105 mm, e.g. 15025, 15026, 15029, 15030, 15032, 15033, polypropylene, suitable for 13034		

Part No.	Description	Max. speed (rpm)	Max. gravitational field (x g)
14035	Round carrier for 1 tube with flat bottom or skirt approx. 30 ml, max. Ø 25/31 x 65 – 95 mm incl. cap, e.g. sterilin tubes, graduated up to 20 ml, see www.bibby-sterilin.co.uk , no. 03008, polypropylene, suitable for 13034		
14032	Round carrier for 1 round bottom tube 50 ml, max. Ø 35/38 x 70 - 105 mm, e.g. 15049, 15050, 15056, polypropylene, suitable for 13034		
14033	Round carrier for 4 round bottom tubes, max. Ø 13.5/17.5 x 70 - 110 mm, e.g. Vacutainer, polypropylene, suitable for 13034		
14034	Round carrier for 3 round bottom tubes 10-15 ml, max. Ø 17.3/19 x 80 - 110 mm, e.g. 15015, 15020, 15022, 15023, 15024 and Monovettes 9 and 10 ml, polyethylene, suitable for 13034		
13031	Bucket, aluminium, incl. polysulfone screw cap 17151, for 1 culture tube 50 ml (Nunc, Corning, Falcon, Greiner) e.g. 15151, suitable for 11030, 11031, 11032, 11035, 11036, max. radius 14.4 cm, min. radius 4.2 cm	4 000	2 576
13060	Adapter for 1 culture tube 15 ml (Nunc, Corning, Falcon, Greiner) e.g. 15115, suitable for 13031		
13038	Carrier, aluminium, for 2 culture tubes 15 ml (Nunc, Corning, Falcon, Greiner) e.g. 15115, suitable for 11030, 11031, 11032, 11035, 11036, max. radius 14.3 cm, min. radius 5.1 cm	4 000	2 558

Part No.	Description	Max. speed (rpm)	Max. gravitational field (x g)
11121	Swing-out rotor for microtitre plates, incl. 2 buckets 13223, radius edge 11.9 cm, radius max. 10 cm, radius min. 6.5 cm, max. plate height 50 mm	4 000	2 129 1 789 1 163
12011	Angle rotor 6 x 15 ml incl. buckets 13011 for round or conical tubes, max. Ø 17 x 130 mm, e.g. 15015, 15020, 15023, 15024 and culture tubes 15 ml (Nunc, Corning, Falcon, Greiner) e.g. 15115 and for 6 reaction vials 1.5-2.2 ml, e.g. 15008, 15040, max. radius 9.6 cm, min. radius 3 cm (with 100 mm tube length), angle 35°	4 000	1 717
12012	Angle rotor 8 x 15 ml incl. buckets 13011 for round or conical tubes, max. Ø 17 x 115 mm, e.g. 15015, 15020, 15023, 15024 and for 8 reaction vials 1.5-2.2 ml, e.g. 15008, 15040, max. radius 9.6 cm, min. radius 3 cm (with 100 mm tube length), angle 35°	4 000	1 717
12061	Angle rotor 30 x 15 ml incl. buckets 13011 for round or conical tubes, max. Ø 17 x 120 mm, e.g. 15015, 15020, 15023, 15024, Monovettes and culture tubes 15 ml (Nunc, Corning, Falcon, Greiner) e.g. 15115, 2 lines, angle 33°; max. radius 13.9 cm, min. radius 7.8 cm, max. radius 11.6 cm, min. radius 5.8 cm	4 000 4 000	2 486 2 075
12062	Angle rotor 20 x 15 ml incl. buckets 13011 for round or conical tubes, max. Ø 17 x 120 mm, e.g. 15015, 15020, 15023, 15024, Monovettes and culture tubes 15 ml (Nunc, Corning, Falcon, Greiner) e.g. 15115, max. radius 13.9 cm, min. radius 7.8 cm, angle 33°	4 000	2 486

Part No. Description

Adaptors, plastic vessels and glass tubes

13021	Adapter for PCR-tube 0.2 ml, Ø 5.85/6.95 x 20/23.4 mm, suitable for 12011, 12012, polyallomer
13000	Adapter for reaction vials 0,25-0,4 ml 15014, suitable for 12011, 12012, polyallomer
13002	Adapter for reaction vials 0.5-0.75 ml 15005, Ø 7.9/10 x 28/31 mm, suitable for 12011, 12012, polyallomer
15014	Reaction vials 0.4 ml (Beckman system), polypropylene, 1 pack contains 100 pcs., suitable for 13000
15005	Reaction vials 0.5 ml, Ø 7.9/10 x 28/31 mm, 1 pack contains 100 pcs., suitable for 13002
15008	Reaction vials 1.5 ml, 1 pack contains 100 pcs., suitable for 12011, 12012
15040	Reaction vials 2.2 ml, 1 pack contains 100 pcs., suitable for 12011, 12012
15060	Polystyrene tube 5 ml (RIA tube), Ø 12 x 75 mm, suitable for 13033, 14033
15007	Centrifuge glass tube 7 ml, Ø 12 x 100 mm, suitable for 14029, 17907
15027	ditto, graduated, 0 - 5.5 ml in steps of 0.1 ml
15015	Centrifuge glass tube 10-15 ml, Ø 16 x 100 mm, suitable for 12011, 12012, 12061, 12062, 13020, 13035, 14030, 17215
15024	ditto, graduated, 0 - 10 ml in steps of 0.1 ml
15020	Polystyrene tube 15 ml, Ø 17 x 100 mm, suitable for 11031, 11035, 12011, 12012, 12061, 12062, 13020, 13035, 14034, 17215
15021	Polypropylene stopper for 15020, 15023

Part No.	Description
15023	Polypropylene tube 15 ml, Ø 17 x 100 mm, suitable for 12011, 12012, 12061, 12062, 13020, 13035, 14034, 17215
15025	Centrifuge glass tube 25 ml, Ø 24 x 100 mm, suitable for 14031, 17225, 17925
15026	ditto, graduated, 5 - 25 ml in steps of 1 ml
15050	Centrifuge glass tube 50 ml, Ø 34 x 100 mm, suitable for 14032, 17250, 17950
15056	ditto, graduated, 4 - 50 ml in steps of 1 ml
15049	Polycarbonate tube 50 ml, Ø 34 x 100 mm, graduated 0 - 50 ml in steps of 1 ml, suitable for 14032, 17250, 17950
15100	Centrifuge glass tube 100 ml, Ø 44 x 100 mm, suitable for 11032, 13099
15106	ditto, graduated, 1 - 100 ml in steps of 1 ml
15102	Polypropylene tube 100 ml, Ø 45 x 100 mm, suitable for 13099
15103	Polycarbonate tube 100 ml, Ø 45 x 100 mm, graduated 2 - 100 ml in steps of 2 ml, suitable for 13099
15115	Culture tube with screw cap 15 ml, pointed bottom, suitable for 12061, 12062, 13060, polypropylene
15151	Culture tube with screw cap 50 ml, pointed bottom, suitable for 13031, polypropylene



Spare parts

13011	Nylon bucket for 1 tube 15 ml, suitable for 12011, 12012, 12061, 12062
17130	Round sealing cap, polysulfone, clear for 13034
17151	Round sealing cap, polysulfone, clear for 13031

Further accessories available on request.

Maximum speed for tubes

Some tubes, e.g. centrifuges glass tubes, microtubes, culture tubes, polyflor tubes and especially high volume tubes can be used in our rotors, buckets and adapters at higher speeds than their breaking limit. We recommend to always fill up the tubes and to follow the recommendations of the manufacturer.



1.3 Scope of Supply

The following belongs to the centrifuge:

Connection cable	Part No. 269 010
Rotor wrench	Part No. 930 050
20 ml slushing oil	Part No. 70104
Grease for load bearing bolts	Part No. 70284
Spare fuses	Part No. 70149 for 230 V
	Part No. 70102 for 100-120 V

Documentation:

Operating Manual
"Rotor and Accessories, Operation and Use"
EU-Statement of Conformity
Equipment Decontamination Certificate

Accessories according to your order, our order confirmation and our delivery note.

Rotor Part No.	Rotor No.
.....
.....
.....
.....
.....
.....
.....

1.4 Standards and Regulations

Please refer to the enclosed EU-Statement of Conformity.











1.5 Safety Instructions

According to the German trade association regulation BGR500 chapter 2.11 part 3 the owner of the instrument is advised to take care of the following points:

1. According to BGR500 the owner has to provide operating instructions based on those of the manufacturer and to inform the employees accordingly.
2. For safety reasons these operating instructions must clearly state that the stamped max. speed of the used rotor and the max. allowable filling quantity must not be exceeded.
3. If the density of the material exceeds 1.2 g/cm^3 , the max. speed of the centrifuge must be reduced (see formula chapter 7.1.2).
4. Operation of the centrifuge in hazardous locations is not allowed.
5. During operation the centrifuge must not be moved. Leaning against or resting on the centrifuge is not allowed.
6. Do not spin explosive or highly inflammable materials.
7. Substances which could damage the material of the centrifuge, the rotors or the buckets anyhow must not be centrifuged or only under consideration of special safety measures. Infectious, toxic, pathogene or radioactive substances must be centrifuged in certified rotors only.
8. Keep a clearance of at least 30 cm around the centrifuge. Dangerous materials of any kind must not be put down or stored in that area.
9. Attention!
Defective lid relieving devices could cause the centrifuge lid to fall down (contact Service). Risk of bruising!

1.6 Symbol Table

International symbols used for the centrifuge:

Symbol	Title
	Gefährliche elektrische Spannung Dangerous voltage Courant haute tension
	Achtung, Bedienungsanleitung beachten Attention, consult accompanying documents Attention, consulter les documents joints
	Ein (Netzverbindung) On (Power) Marche (mise sous tension)
	Aus (Netzverbindung) Off (Power) Arrêt (mise hors tension)
	Schutzleiteranschluß Protective earth (ground) Liaison à la terre
	Erde Earth (ground) Terre
	Netzstecker ziehen Unplug mains plug Tirer la fiche de prise
	Vorsicht Quetschgefahr Caution! Risk of bruising Attention! Danger de blessure
	Drehrichtungspfeil Arrow direction of rotation Flèche sens de rotation
	Heiße Oberfläche Hot surface Surface chaude

2.1 General Outlay

The new generation of SIGMA laboratory centrifuges is equipped with newest state-of-the-art electronics and is driven by brushless, silent and long-life asynchronous motors.

The problem of carbon brush change is no longer existent and as there is no carbon dust pollution, operation in clean rooms is possible if the appropriate accessories are used.

2.2 Construction and Constructive Safety Measures

The centrifuge is built into a solid steel housing. The centrifuge lid is also made of solid steel. From the back, the lid is secured by solid hinges and at the front by two separate cover locks.

2.3 Drive

The drive motor is a well dimensioned asynchronous motor.

2.4 Operation and Display

The graphical LCD display is hermetically sealed, operation is executed via two knobs. Any operating status is indicated.

2.5 Electronics

The electronics controlled by a microprocessor allows extensive adaptations of the centrifuge to the different tasks. The following parameters can be set:

- Speed in steps of 100 rpm
- RCF in steps of 250 x g
- Time in steps of 1 min (max. 30 min)
- Continuous operation
- Short-time operation

2.6 Safety Devices

Apart from the passive safety devices due to the instrument's mechanical design there are the following active precautions for your safety:

2.6.1 Lid Lock, Cover Closing Device

The centrifuge can only be started when the lid is correctly closed. The cover locks must close. The lid can only be opened when the rotor has stopped. If the lid is opened by the emergency release during operation, the centrifuge will immediately switch off and decelerate brakeless. If the lid is open, the drive is completely separated from the mains supply, that means starting of the centrifuge is impossible (refer to point 7.2.4 "Emergency lid release").

2.6.2 Standstill Monitoring

Opening of the centrifuge lid may only be possible, if the rotor is at standstill. This standstill is checked by the microprocessor.

2.6.3 System Check

An internal system check monitors data transmission and the sensor signals with regard to plausibility. In case of a malfunction an error message is displayed under a number in the time area.

2.6.4 Ground Wire Check

For ground wire check there is a ground screw at the rear panel of the centrifuge. A ground wire check can be carried out using an appropriate measuring instrument.

2.6.5 Imbalance Monitoring System

The display "Imbalance" indicates by lightning up that the centrifuge is in the still permitted imbalance area.

In the event that bigger uneven loading leads to imbalance, the drive is switched off during acceleration or during run.

"Imbalance" and the sequence of bars are flashing. The lid opens automatically after the rotor has come to a standstill. In both cases the loading has to be checked and balanced.

3.1 Unpacking of the Centrifuge

Open cardboard. Take out the box containing accessories. Remove upper foam cushions. Lift centrifuge upwards with a lifting device or with several persons. When lifting or carrying the centrifuge please always reach under the instrument from the side.

Attention: The instrument is heavy!

Please keep case for possible transport of centrifuge later.

3.1.1 Transport Safety Device

The SIGMA 2-5 has a transport safety device which must be removed before start-up.

Procedure:

Remove the foamed plastic ring from the centrifuge chamber.

The transport safety device should be kept for possible transport of the centrifuge (service, repair).

3.2 Installation

3.2.1 Site

All energy consumed by the centrifuge is converted into heat and emitted into the ambient air. Therefore, sufficient ventilation is important. As the air-ducts in the unit must be open, keep a clearance of at least 30 cm around the centrifuge. Also, the centrifuge shouldn't be positioned near radiators and should not be directly exposed to sunshine.

The table should have a solid, even top.

For normal operation the ambient temperature should not fall below 10 °C and not exceed 35 °C. The max. humidity of air is 80 %. During transport from cold to warmer places water will condensate inside the centrifuge. It is important that there is enough time for drying before the centrifuge can be started again.

3.2.2 Connection

The operating voltage on the name plate must correspond to the local supply voltage!

SIGMA laboratory centrifuges are units of safety class I, DIN VDE 0700, and include a three wire power cord 2,5 m long with shockproof right angle plug.

At the back, next to the mains supply, there is an additional ground wire connection where a separate ground wire can firmly be connected to ground any non-hazardous leakage current. The leakage current is harmless but secondary effects would occur.

3.2.3 Fuses / Emergency Circuit Breaker on Site

The centrifuges must be protected typically with 16 A time lag low fuses.

An emergency circuit breaker to cut the power to the centrifuge in the event of a malfunction is required on site. This switch should be located away from the centrifuge, preferably outside the room where the centrifuge is used or at the exit of this room.

3.3 Installation of Rotors and Accessories




1. Open centrifuge lid by pressing the Lid-key.
2. Unscrew rotor tie-down screw from motor shaft (anticlockwise).
3. Lower the rotor straight down onto the motor shaft.
4. Tighten the tie-down screw (clockwise) with the rotor wrench with approx. 5 Nm.

In the event of frequent use the tie-down screw must be loosened by some turns and fastened again. **This should be done once a day or after approx. 20 cycles.** This ensures a proper connection between rotor and shaft (please refer to chapter 6.2 "Care and cleaning of accessories" as well).

5. Use only appropriate vessels for the rotor (please refer to chapter 1.2 "Suitable accessories" as well).
6. Fill vessels external to the centrifuge.
7. Put or screw on covers of vessels.
8. Opposite places of the rotors must always be loaded with same accessories and same filling.
9. In angle rotors the plastic vessels must always be totally filled to avoid cracks of vessels and leakages or loosening of the caps in case of partial filling.

Attention, follow the special comments of chapter 1.5.

10. **Attention:** The centrifuge will absorb smaller differences in weight when loading the rotors. But it is recommended to balance the vessels as accurately as possible in order to ensure a run with minimal vibrations.
11. Don't fix the rotor screw without a rotor. Otherwise you'll destroy the shaft.
12. Rotors with lid should always be run with their lid. The rotor lid is tightened by the lid screw or by hand if you use a lid with a knurled screw. Correct fastening must be ensured. **Attention: The lid screw serves for fastening of the lid onto the rotor only, not for fastening of the rotor onto the drive!** Before installation of the lid, the correct fastening of the rotor fixing screw must always be checked using a wrench.



3.4 Initial Start-Up

Attention!

Before initial start-up please take care that your centrifuge is orderly installed (refer to chapter 3.2 "Installation").

3.4.1 Switching on of the Centrifuge

Press mains switch (right side panel).

- The centrifuge display is illuminated.

3.4.2 Opening Lid




Press the Lid-key.

- The lid opens.

3.4.3 Installation of an Angle Rotor and a Swing-out Rotor

Put a rotor onto the shaft and fasten it by screwing the rotor tie-down screw clockwise onto the drive shaft. Please use the supplied rotor wrench and hold the rotor with one hand (refer to chapter 3.3 "Insertion of rotor and accessories").

4.1 Operating Panel



The centrifuge is operated via the operating panel. Keys can be pressed when their LED is on.

When power is applied, all LEDs and the display are illuminated for a short time.

4.1.1 Start-key



This key can be used for the following:

- starting centrifuge operation,
- terminating a previously started deceleration process and restarting the centrifuge,
- shifting to short-run (refer to chapters 4.2.3.2.1, 4.2.3.2.2).
- saving a program. The parameters speed and time preselection can be locked, thereafter a change is impossible. Starting and stopping the centrifuge and opening of the lid is always possible.
Open centrifuge lid. Press the start-key 3 x and hold for approx. 2 seconds when pressing for the third time. A flashing bar appears in the time and/or speed panel. Cancellation of this function occurs in the same way.

The centrifuge can be started when

- the lid is closed,
- the Start-key is illuminated.

4.1.2 Stop-key with Soft and Softstop function




This key can be used

- to early terminate the run,
- to activate the softstop function. This means a doubling of the deceleration time.
- to activate the soft function resulting in double the braking and starting time.

An individual function softstart is not possible.

This is possible when

- the centrifuge stands still:
Activate or deactivate "soft stop" or "soft" by pressing the Stop-key.
- the centrifuge runs, Stop-key is illuminated:
Press Stop-key. During the deceleration process "soft stop" or "soft" can be activated or deactivated by pressing the Stop-key again. In case of activation "soft stop" or "soft" is indicated in the display.



The centrifuge decelerates to a complete stop. Deceleration can be terminated by pressing the Start-key again and the centrifuge can be restarted.

4.1.3 Lid-key



This key is used to open the lid.

This can only be executed if

- the mains switch is "ON",
- the Lid-key is illuminated.

A flashing Lid-key indicates that only on lid lock has locked. Press down left and right side of the lid to lock both lid locks.


The lid opens automatically after standstill of the rotor!

4.1.4 Knobs

The parameters are selected by turning the knobs.

The time values are changed by turning the left knob, the right knob is used for changing the speed/RCF values.

4.2 Displays



The display consists of two areas:

The left area is the time display, the right one the speed/RCF display.


4.2.1 Speed/RCF (Relative Centrifugal Force)

Rotor dimensions and speed determine the RCF value. Entry of one value automatically determines the other one.


4.2.1.1 Speed

The upper scale of the right display area indicates rpm x 1000. This value can be set by turning the right knob.

After entry of a speed a single bar indicates the set speed.




During the run the actual speed is indicated by a sequence of bars.



The speed value can be changed during the run by turning the right knob.

If during the run the speed is set below the actual speed, the new set speed is indicated by a flashing bar and the drive is decelerating to the new value.




4.2.1.2 Relative Centrifugal Force (RCF)


The relative centrifugal force (RCF) is the acceleration which the sample is exposed to.

The lower scale of the right area indicates the RCF value x 1000. This can be set by turning the right knob.

After entry of a RCF value a single bar indicates the set speed.




A sequence of bars indicates the actual RCF value during a run.



The RCF value can be changed during a run by turning the right knob.

If during the run the RCF value is set below the actual value, the new set RCF value is indicated by a flashing bar.



Please refer to chapter 1.2 "Suitable accessories" for the maximum values for the used rotor combination. For this example rotor combinations 11035 and 11032 have been chosen.

4.2.2 Time

Depending on the run mode this area indicates the total set time, the remaining run time or the elapsed time.


The time value as well as different run modes can be set by turning the left knob, even during centrifugation.

If the time is changed during centrifugation, the total newly entered period is run. the time elapsed before is not considered.


4.2.3 Standard Run Mode

For the standard run mode the time is set in minutes by turning the left knob ("min" is displayed).


This total run time (max. 30 minutes) is displayed by a sequence of bars before starting the centrifuge.



During centrifugation this time is counted down, the sequence of bars indicates the remaining run time. The set time is indicated by one single bar.




The centrifuge run can be early terminated by pressing the Stop-key. The set time and the remaining run time remain on the display.




4.2.3.1 Continuous Run Mode

During continuous run the centrifuge accelerates up to the set speed and the run must be terminated manually.

Select ∞ turning the left knob clockwise and press the Start-key to activate the continuous run.




Unlike the standard run mode the minutes are counted up. A sequence of bars indicates the elapsed run time. One single bar at ∞ is indicating the continuous run mode.



After 30 minutes the remaining run time is no longer displayed but the run continues.

The continuous run is terminated by pressing the Stop-key and the centrifuge decelerates to a complete stop. The elapsed time is displayed in minutes.



The continuous run mode can be left by selecting a set value as well.

4.2.3.2 Short-Time Operation


There are two options:

4.2.3.2.1 Short-Run 1


Press the Start-key for the whole short-run. The centrifuge accelerates with maximum power to the maximum speed (3900 rpm). After release of the Start-key the centrifuge decelerates with maximum power to standstill.

During short-run the time is counted up in seconds, "sec" is displayed.

A sequence of bars is indicating the elapsed run time. One single bar at "short run" indicates that short-time operation had been selected. A flashing bar indicates the maximum speed.




After termination of the short-run the elapsed run time in seconds remains on the display.




4.2.3.2.2 Short-Run 2

During this short-run the centrifuge is running max. 30 seconds and accelerates to the set speed.

Select "short run" by turning the left knob and press the Start-key.




During the short-run the time is counted up in seconds. A sequence of bars indicates the elapsed run time. One single bar at "short run" indicates that short-time operation had been selected.



After 30 seconds the centrifuge decelerates to a complete stop.

The short-run can be early terminated by pressing the Stop-key. The centrifuge decelerates to standstill. The elapsed run time in seconds remains on the display.




The "short run" can be left by selecting a set value as well.

5.1 Practical Notes for Centrifugation


1. Locate centrifuge horizontally on a level surface.
2. Ensure safe location.
3. Keep at least 30 cm free space around the centrifuge.
4. Provide for sufficient ventilation.
5. Tighten rotor firmly onto motor shaft.
6. Avoid imbalance.
7. Load opposite places with same accessories.
8. Centrifugation with low capacity:

An example is the fixed angle rotor 24 x 2,2 ml and the rotor 8 x 15 ml which should be loaded with 4 vessels only.

The vessels should be placed symmetrically so that the rotor is loaded evenly. Loading e.g. only one position is not allowed.



9. Load all positions of swing-out rotors.
10. Load vessels outside the centrifuge.
11. Fill vessels carefully to same weight. Imbalances would result in increased wear of bearings.
12. Use perfect accessories only.
13. Avoid corrosion to accessories by careful maintenance.
14. Spin infectious material in sealed rotors and buckets only.
15. Do not spin explosive or highly inflammable materials.
16. When centrifuging substances with a density $> 1,2 \text{ g/cm}^3$ the allowable max. speed must be reduced (refer to chapter 7.1.2 "Density").
17. Grease the bolts of the rotor and the nutes of the buckets slightly with grease for load bearing bolts.



5.2 Forbidden Centrifuging Operations

1. Operation of not carefully installed centrifuge.
2. Operation without front or back panels.
3. Operation by non authorized personnel.
4. Operation with rotor not installed properly (refer to chapter 3.3).
5. Operation with overloaded rotors.

The load for a rotor is limited by its design and the max. speed (see rotor/bucket engraving) and must not be exceeded. The rotors are intended for liquids of max. homogeneous density of 1.2 g/cm^3 if centrifuged at max. speed. If liquids of higher density are used, the speed must be reduced accordingly (refer to chapter 7.1 "Mathematical relations").

6. Operation with rotors and adapters showing corrosion or other defects.
7. Operation of very corrosive substances which can cause damages to material and affect mechanical strength of rotors and adapters.
8. Operation of rotors and accessories not allowed by the manufacturer. The use of poor commodity goods is not recommended. At high speeds breaking glass or bursting vessels can cause dangerous imbalances.
9. Operation in hazardous locations.
10. Operation with vessels of improper size.
11. Centrifugation of improper material.
12. Operation with partially filled plastic tubes in high-speed angle rotors.
13. Lifting or moving of the centrifuge during operation. Leaning against or resting on the centrifuge is not allowed.
14. Do not place potential dangerous material - e.g. glass vessels containing liquids - near the centrifuge.
15. Attention: Do not open cover and/or reach into rotor chamber unless the rotor is at standstill. Never attempt to override the lid interlock system while the rotor is spinning.
16. Such materials are prohibited which chemically interact vigorously.
17. Do not spin explosive or inflammable materials.

18. Substances which could damage the material of the centrifuge, the rotors or the adapters must not be centrifuged. Infectious, toxic, pathogene or radioactive substances must be centrifuged in certified rotors and vessels only and all necessary safety precautions are taken.
19. Operation with incompletely loaded swing-out rotor or angle rotor with interchangeable buckets.

A rotor must always be loaded completely, empty places are not allowed! Opposite buckets or carriers may nevertheless be empty. Mixed loading is allowed, if opposite places are loaded with same buckets and carriers of same weight.

6.1 Care and Cleaning of Centrifuge

Please use water-soluble, mild detergents for cleaning. Avoid corroding and aggressive substances. Do not use alkaline solutions or solvents or agents with abrasive particles. Before using detergents or decontamination agents which had not been recommended by us, the user has to contact us to make sure that such procedure would not damage the centrifuge.

Remove product particles from the rotor chamber using a cloth or paper towel. It is recommended to open the cover when the centrifuge is not in use so that moisture can evaporate. Increased wear of the motor bearings will thus be avoided. **If there is the risk of toxic, radioactive or pathogene contamination, special safety measures must be kept.**

6.2 Care and Cleaning of Accessories

For care of accessories special safety measures must be considered as these are measures ensuring operational safety at the same time.

Chemical reactions as well as stress-corrosion (combination of oscillating pressure and chemical reaction) can affect or destroy the metals. Hardly detectable cracks on the surface expand and weaken the material without visible signs. When detecting a visible damage of the surface, a crack, a mark or any other change, also corrosion, the part (rotor, etc.) must be replaced immediately.

In order to avoid corrosion, rotor incl. tie-down screw and cover seal and adapters must be cleaned and greased regularly with the supplied slushing oil (Sigma part no.: 70104 for 20 ml slushing oil). Before using detergents or decontamination agents which had not been recommended by us, the user has to contact us to make sure that such procedure would not damage the centrifuge. The rotor tie-down screw must be greased using grease (Sigma part no.: 70284).

Cleaning of accessories should be done outside of the centrifuge once a week or preferably after every use. Adapters should be removed. After this the parts should be dried with a soft cloth or, alternatively, in a drying chamber at approx. 50 °C. **If there is the risk of toxic, radioactive or pathogene contamination, special safety measures must be kept.**

Especially aluminium parts are extremely corrosive. A neutral cleaning detergent with a pH-value between 6 and 8 should be used for such parts. Alkaline agents exceeding pH 8 must be avoided. Especially aluminium parts must be greased regularly with slushing oil. This procedure essentially increases life time and reduces corrosion.

Careful maintenance increases life time and avoids premature failure of the rotor. Corrosion or resultant damages which are caused by insufficient care do not constitute a warranty claim.

6.2.1 Service notes

The rotor part no. 12034 consists of polypropylene. The life essentially depends on compliance with the chemical resistance tables of the producer (e.g. Bayer, BSAF).

The centrifuge and the rotors should not be exposed to intensive UV radiation and longer thermal stress. Cleaning should be done using mild detergents.

If required the rotor can be removed. After reinstallation please tighten the rotor tie-down screw properly.

In the event of material changes (cracks) or deformation and uneven run the rotor must no longer be used and the manufacturer should be informed accordingly.

With increasing temperature the chemical resistance of plastics will reduce.

The chemical resistance must be checked beforehand!

Careful maintenance increases life time and avoids premature failure of the rotor. Corrosion or resultant damages which are caused by insufficient care do not constitute a warranty claim.

6.3 Rotor Trunnion Pins

The trunnion pins of the rotor should always be greased as only this ensures evenly swinging of buckets and thus quiet run of the centrifuge (part no. 70284 Grease).

6.4 Glass Breakage

In case of glass breakage all glass particles must be carefully removed. Rubber inserts have to be cleaned carefully and possibly be replaced. If a problem has occurred, the following has to be considered:

Glass particles in the rubber cushion will cause glass breakage again.

Glass particles in the rotor pins prevent buckets and carriers from swinging evenly which will cause an imbalance.

Glass particles in the centrifuge chamber will cause metal abrasion due to the strong air circulation. This dust will not only pollute the centrifuge chamber, the rotor and the material to be centrifuged but also damage the surfaces of the accessories, the rotors and the centrifuge chamber.

In order to totally remove the glass particles and the metal dust from the rotor chamber, it is advisable to grease the upper part of the centrifuge chamber with e.g. Vaseline. Then the rotor should rotate for some minutes at a moderate speed. The glass and metal particles will now collect at the greased part and can easily be removed with a cloth together with the grease. If necessary repeat this procedure.

6.5 Sterilization and Disinfection of Rotor Chamber and Accessories

All usual disinfectants like eg. Sagrotan, Buraton or Terralin (to obtain at chemist's shops) can be used. The centrifuges and the accessories consist of different materials. A possible incompatibility must be considered. Before using detergents or decontamination agents which had not been recommended by us, the user has to contact us to make sure that such procedure would not damage the centrifuge. For sterilization by steam resistance to temperature of the individual material must be checked (refer to point 6.5.1 "Autoclaving"). Please contact your laboratory safety officer regarding proper methods to use. **If dangerous materials are used, the centrifuge and the accessories must be disinfected.**

Principally we want to point out that for centrifuging of e.g. infectious material certified and hermetically sealed accessories have to be used in order to avoid that the centrifuge is contaminated.

6.5.1 Autoclaving

The life of the accessories essentially depends on the frequency of autoclaving and use. When the parts are showing changes in colour or structure or in the event of leaks etc., the accessories have to be replaced.

During autoclaving the caps of the tubes must not be screwed on to avoid deformation of the tubes. It can not be excluded that plastic parts, e.g. lids or carriers, would deform during autoclaving.

Autoclaving:

Accessories	max. temp. °C	min. time min	max. time min	max. cycles
Glass tubes	134-138	3	5	-
Polycarbonate tubes	115-118	30	40	20
Polypropylene tubes	115-118	30	40	30
Teflon tubes	134-138	3	5	100
Aluminium rotors	134-138	3	3	-
Polypropylene rotor 12034	115-118	30	40	20
Polypropylene rotor 12124	115-118	30	40	20
Polycarbonate/Polyallomer lids for angle rotors	115-118	30	40	20
Polysulfone lids for angle rotors	134-138	3	5	100
Aluminium buckets	134-138	3	5	-
Polycarbonate caps for buckets	115-118	30	40	50
Polypropylene caps for buckets	115-118	30	40	50
Polysulfone caps for buckets	134-138	3	5	100
Rubber adapters	115-118	30	40	-
Rubber cushions	115-118	30	40	-
Round carriers for 13104/ 13117, Polypropylene	115-118	30	40	-
ditto, Polyallomer and Polycarbonate	115-118	30	40	-
Round carriers for 13350/ 13550, Polypropylene	115-118	30	40	-
Rectangular carriers, Polypropylene	115-118	30	40	-
ditto, Polyallomer and Polycarbonate	115-118	30	40	-

6.6 Checks by Operator

The operator has to ensure that no important part of the centrifuge is damaged. This especially refers to:

1. Motor suspension
2. Concentricity of the motor shaft
3. Fastening of the trunnions in the rotor
4. Rotors and accessories have a limited life. For safety reasons a regular check is recommended after 50.000 cycles. Any changes like e.g. corrosion, cracks, material abrasion etc. require special attention.
5. Screw connections have to be tight.

Furthermore, the earth wire must be checked regularly.

7.1 Mathematical Relations

7.1.1 Relative Centrifugal Force (RCF)

The parameters speed, RCF and the diameter of the rotor are interrelated via the following formula:

$$\text{RCF} = 11,18 \times 10^{-6} \times r \times n^2$$

If two values are given, the third value is determined by the equation. If the speed or the rotation radius is changed, the resulting RCF will be recalculated. If the RCF is altered, the speed under consideration of the radius is adapted accordingly.

r = radius in cm
n = speed in rpm
RCF without dimension

7.1.2 Density

If the density of the liquid is higher than 1.2 g/cm³, the allowed maximum speed of the centrifuge is calculated according to the following formula:

$$n = n_{\text{max}} \times \sqrt{(1,2 / \text{Rho})}$$

Rho = density in g/cm³

7.2 Error Correction

Most of the errors can be reset by power off/on. In the event of a short power failure during a run, this run is interrupted and can be continued by pressing the Start-key.

No indication on the display: Actions:

- | | |
|--|--|
| <ul style="list-style-type: none">– No power in the socket?– Power cord plugged in and line voltage present?– Input fuse ok?
– Power switch on?– Lid closed? | <ul style="list-style-type: none">• Check fuse in mains supply.• Plug in power cord correctly.
• Replace input fuse (see nameplate for rating).• Switch on power.• Close lid (please refer to chapter 7.2.3 "Lid cannot be opened/closed."). |
|--|--|

7.2.1 Centrifuge cannot be Started

- | | |
|--|---|
| <ul style="list-style-type: none">– Start-key LED not illuminated:

– Lid-key LED flashing: | <ul style="list-style-type: none">• Power off/on. If error occurs again, call service.

• Open and close the lid again. If error occurs again although both locks engaged, call service. |
|--|---|

7.2.2 Centrifuge Decelerates during Operation

- | | |
|--|--|
| <ul style="list-style-type: none">– Centrifuge displays an error 1 to 11 after power on. | <ul style="list-style-type: none">• Power off/on. If error occurs again, call service (please refer to chapter 7.3.1 "Error codes"). |
|--|--|

7.2.3 Lid cannot be Opened/Closed

- When first trying to open the lid the locks are not released. The Lid-key LED is flashing. Open and close lid again.

- Check/clean the lid seal. Apply talcum powder to the seal to avoid renewed sticking.

7.2.4 Emergency Lid Release

In the event of e.g. a power supply failure it is possible to manually open the lid. At the bottom panel there are two stoppers which can be removed e.g. with a screw driver. Before, unscrew the screw that is in the stopper. The lid can be released by pulling the visible strings. Thereafter, stoppers and screws must be fixed again.

Attention!

The lid may only be unlocked and opened when the rotor is at standstill.

7.2.5 Problems with the Centrifuge?

Please contact your supplier for support or in the event of malfunctions and for supply of spare parts.

7.3 Error Mode

In the error modes the Start-, Stop- and Lid-keys and the bars in the speed display are flashing. The error codes are displayed on the time display.

7.3.1 Error Codes

Error no.	Kind of error	Actions
1	Tacho signal disturbed	• Power off/on
4	Lid does not open after pressing the Lid-key three times	• Power off/on • Emergency lid release
2, 3, 5 - 11	Internal fault	• Power off/on

Should it not be possible to repair the failure, please contact Service!

7.4 Speed-RCF-Diagram

An additional help is the enclosed Speed-RCF-Diagram.

7.5 Declaration of Decontamination / Return Declaration

The following declarations serve for keeping safety and health of our employees. Fill in the forms and attach them when returning centrifuges, accessories and spare parts. Please understand that we cannot carry out any work before we have the declarations. (We recommend to make **several copies of this page.**)



!!! Attention – This form must be glued on outside of the box !!!

Return declaration

	YES	NO
Decontamination declaration inside :		
Unit / component contaminated :		
Unit / component unused (new) :		

!!! Attention – This form must be glued on outside of the box !!!



Please make some copies before removing this page.

✂-----
Declaration of Decontamination of Centrifuges, Accessories and Spare Parts

This declaration may only be filled in and signed by authorised staff.



Repair Order dtd. : _____

Order No. : _____

Type of unit : _____ Serial No. : _____

Type of unit : _____ Serial No. : _____

Type of unit : _____ Serial No. : _____

Type of unit : _____ Serial No. : _____

Accessories : _____

Is the equipment free from harmful substances ? YES NO

If not, which substances have come into contact with the equipment?

Name of the substances : _____

Remarks (e.g to be touched with gloves only) ; _____

General characteristics of the substances :

Corrosive	<input type="radio"/>	Explosive	<input type="radio"/>
Biologically hazardous	<input type="radio"/>	Radioactive	<input type="radio"/>
Toxic	<input type="radio"/>		

In combination with which substances may hazardous mixtures develop?

Name of the substances : _____

Has the equipment been cleaned before shipment? YES NO

Is the equipment decontaminated and not harmful to health? YES NO

Prior to repair, radioactively contaminated components must be decontaminated according to the valid regulations for radiation protection.

Legally Binding Declaration

I / we hereby declare that the information on this declaration are correct and complete.

Company / Institute : _____

Street : _____

Postcode, City : _____

Tel. : _____ FAX : _____

Name : _____

Date : _____ Stamp : _____

Signature : _____

✂-----

Please make some copies before removing this page.

